
User's manual
Manuel d'utilisation
Manual de uso Benutzerhandbuch
Manuale d'uso
Handleiding voor de gebruiker
Manual de utilização
Användningshandbok
Εγχειρίδιο χρήσης

PIEZOTOME
Piezo.Ultrasonic.Surgery.Unit
SOLO

PIEZOTOME
Piezo.Ultrasonic.Surgery.Unit
SOLO LED

SATELEC
ACTEON

ENGLISH	4
FRANÇAIS	20
ESPAÑOL	36
DEUTSCH	52
ITALIANO	68
NEDERLANDS	84
PORTUGUÊS	100
SVENSKA	116
ΕΛΛΗΝΙΚΑ	132

1	INTRODUCTION	5	11	SERVING AND STERILIZATION	
				Servicing the device	12
2	WARNINGS	5-6		Servicing the ultrasound handpiece	12
				Servicing tips	12
3	DESCRIPTION			Servicing irrigation lines	12
	Physical description	6-7	12	CHECKING / MAINTENANCE	
	Technical characteristics	7-8		Checking	12
4	INSTALLATION / COMMISSIONING			Maintenance	12
	Unpacking the device	8		Replacing fuses	13
	Recommendations	8		Troubleshooting	14
	Installation	8-9	13	ELECTROMAGNETIC COMPATIBILITY	
	First start-up	9		Electromagnetic emissions	15
5	IRRIGATING SOLUTIONS	9		Electromagnetic immunity	16
				Electromagnetic immunity /	
6	IRRIGATION LINES	9		radiofrequency mobile equipment	17
				Recommended separation distances	18
7	SETTINGS / PARAMETERS / INTERFACE			Cable lengths	18
	Start-up	9	14	DISPOSAL AND RECYCLING	19
	Parameters	9			
	Flush / prime	9	15	MANUFACTURER'S LIABILITY	19
	Irrigation	9			
	Ultrasound function	9-10	16	REGULATIONS	16
	Tool kit	10			
8	SAFETY	10-11	17	SYMBOLS	148-149
9	USING THE DEVICE	11	18	CUSTOMER RELATIONS	
				Manufacturer identification	154
10	SWITCHING OFF THE DEVICE	11		Subsidiaries	154-155

The list of accessories referred to in this manual is not exhaustive. Consult your local dealer for further information.

I - INTRODUCTION

Congratulations on the purchase of your new PIEZOTOME SOLO™/PIEZOTOME SOLO LED device.

To maximise use of this cutting-edge product for as long as possible, please read these instructions carefully before switching on, using or servicing the device.

Sentences containing the symbol  are points of particular importance.

Designed for clinical use:

Created by SATELEC, PIEZOTOME SOLO/PIEZOTOME SOLO LED is a device designed to perform dental surgery including osteotomies, osteoplasties, periodontal and implant surgery...

II - WARNINGS

CAUTION:

United States Federal Law restricts the use of this device solely to qualified, trained and competent dental health practitioners or under their supervision.

The surgical tips designed for PIEZOTOME SOLO/PIEZOTOME SOLO LED are not compatible with the IMPLANT CENTER, and vice versa.

PIEZOTOME SOLO/PIEZOTOME SOLO LED can only take tips adapted to PIEZOTOME 2 or IMPLANT CENTER 2.

To reduce the risk of accidents, the precautions stipulated below must be taken:

Device operators:

- Use of the device is restricted solely to graduate dental health professionals, who are competent and qualified within the usual context of their activities.
- If you have received this device by mistake, please contact the device supplier to arrange for its return.

Interactions:

Interference may occur when used on patients with cardiac pacemakers.

This system emits electromagnetic fields, which means there are some potential risks. The malfunctioning of implantable devices such as cardiac pacemakers and ICDs (implantable cardioverter defibrillator) is possible:

- Ask patients and users if they have an implanted device before using this product. Explain the circumstances to them.
- Weigh the risks and benefits and contact your patient's cardiologist or appropriate qualified healthcare professional prior to performing the treatment.
- Keep this product away from implanted devices.
- Make appropriate emergency provisions and take immediate action if patients become ill.
- Symptom including a raised heartbeat, irregular pulse, and dizziness may signal problems with a cardiac pacemaker or ICD.

Electrical connection:

- Ask an approved dental installation technician to connect your device to the electrical network.
- Warning: To prevent any risk of electric shock, this device shall only be connected to a supply network fitted with a protective ground.
- The electrical network to which the device is connected shall comply with the standards in force in your country.
- If loss of electrical power during use is likely to generate an unacceptable risk, the device must be connected to a suitable power source (e.g. UPS).
- Do not attempt to connect to the device any other accessories than those supplied by SATELEC.

Overheating, burns:

To avoid the risk of burns or of overheating:

- Please refer to the clinical booklet provided with the device.
- The irrigation flow must be used all the time during the clinical act (please refer to the clinical booklet).
- Please replace the tips when the cutting part becomes blunt or when the diamond coated part

becomes smooth and shiny.

- The handpiece must be cleaned properly and regularly as mentioned in the user's manual.

Using the device:

- Do not use the device if it appears to be damaged or faulty.
- Before unplugging the mains power cable, switch the device off.
- To unplug the mains power cable, pull the plug on the end of the cable whilst holding the wall socket with your other hand.
- Never use irrigating solution containers that are not designed to be suspended from the bracket provided.
- The device must only be used with bottles or pouches of sterile water or physiological saline solution.
- The capacity of the irrigating solution container used shall not exceed one litre.
- If you are not using or will be away from the device for a period of time, disconnect the device from the electrical supply network.
- Do not move the device during use.
- Do not apply excessive pressure to the screen.

Environment:

- Do not cover the device or obstruct the air inlets.
- Do not immerse and do not use outside.
- Do not tilt the device at an angle of more than 5°.
- Do not place the device near a heat source.
- Make sure that the cables do not prevent people from moving about freely.
- The device shall be stored in its original packaging, in a suitable place, which does not present a hazard for people.
- The device is not designed to be used in a medical operating room.
- The device is not designed to operate in the presence of anaesthetic gases or any other flammable gas.
- Do not expose the device to water mist or spray.
- Any condensation inside an electrical device is potentially dangerous. - If the device is to be moved from a cold environment to a hot one, it must not be used until it has reached room temperature.
- The device is not designed to operate in the vicinity of ionizing radiation.

- Do not insert metal objects into the device to prevent any risk of electric shock, short-circuit or the emission of hazardous substances.

Servicing:

- Before and after each use, it is vital that you disinfect your device using products recommended by SATELEC.
- Prior to each intervention, it is imperative that you use clean, disinfected and sterilized accessories.

Accessories:

- The device may only be fitted with accessories designed or recommended by SATELEC that are perfectly suited to the purpose.
- Use of accessories from other manufacturers is a potential hazard for you or your patients.

Repair:

- Warning: Do not repair or modify the device without the prior authorization of SATELEC.
- Warning: If the device is modified or repaired, specific inspections and checks shall be performed to ensure that the device can still be used safely.
- In the event of an anomaly, contact the supplier of your device. Do not contact any other repairer who could render your device hazardous for you and your patients.

In the event of doubt, contact an approved dealer or the Customer service department:

- www.acteongroup.com
- E-mail: satelec@acteongroup.com.

III - TECHNICAL DATA

3.1 PHYSICAL DESCRIPTION

The device comprises 5 main parts:

- A console (Fig. 1-1).
- A multifunction pedal (Fig. 4).
- A PIEZOTOME SOLO or PIEZOTOME SOLO LED handpiece and cable (Fig. 1-2).
- A bracket for irrigating solution (Fig. 1-3).
- A mains power cable with ground (Fig. 1-4).

The top of the console incorporates:

- An LCD (Liquid Cristal Display) screen displaying different information (Fig. 1-5).
- A control keyboard enabling you to change the settings of your device (Fig. 1-6).

The front face incorporates the handpiece cable connector (Fig. 1-7).

The right side face incorporates the pump housing (Fig. 1-8).

The back of the box incorporates the following elements:

- 1 mains connector with an earthing pin (Fig. 2-1).
- 1 pedal connector (Fig. 2-2).
- 1 fan (Fig. 2-3).
- 1 bracket support (Fig. 2-4).
- 1 mains isolating switch (Fig. 2-5).
- 1 potential equalizing terminal (Fig.2-6).

3.2 TECHNICAL DESCRIPTION

a) LCD screen / control keyboard

The device can be adjusted by applying moderate pressure to the keyboard buttons.

Button positions (Fig. 3-1 to 5).

Item 1: Increase irrigation flow.

Item 2: Decrease irrigation flow.

Item 3: Flush/prime.

Item 4: Save, validate or select.

Item 5: Programme selection.

The LCD screen displays all information required by the user.

Display positions (Fig. 3-6 to 9).

Item 6: Store in memory.

Item 7: Irrigation flow rate with associated bar graph.

Item 8: Programme with associated bar graph.

Item 9: Device status (ultrasound activated, error detection).

b) Rear face of console

The mains connector (Fig. 2-1) with its earthing pin allows the device to be connected to the electrical network via a removable mains power cable.

The pedal connector (Fig. 2-2) allows the device to be connected to the multifunction control pedal.

The fan (Fig. 2-3) protected by a metal grille ensures the device always operates at optimum performance level.

The bracket support (Fig. 2-4) allows the bracket to be positioned.

The mains isolating switch (Fig. 2-5) is used to switch the device on or off.

A potential equalizing terminal (Fig. 2-6) allows the device to be connected to other installation equipment as required.

c) Right-side face of console

The pump housing (Fig. 1-8) is designed to receive SATELEC irrigation line cassettes.

An irrigation cassette is installed by lifting the pump cover and inserting the cassette horizontally into the purpose-designed housing.

d) Front face

The connector (Fig. 1-7) is designed to receive the handpiece cable connector. Only the PIEZOTOME SOLO or PIEZOTOME SOLO LED handpiece shall be used.

e) Control pedal

The control pedal allows professionals to work in a completely sterile environment by preventing cross contamination risks.

Once the various parameters have been set, the operator no longer needs to use the keyboard.

Pedal definition (Fig. 4-1 to 3).

Item 1: Ultrasound activation.

Item 2: Mode change (from d1 to d4).

Item 3: Arch.

f) Technical characteristics

Manufacturer: SATELEC

Name of device: PIEZOTOME SOLO or PIEZOTOME SOLO LED

Electrical supply:

- Voltage: 100 VAC to 230 VAC

- Frequency: 50 Hz / 60 Hz

- Nominal power: 150 VAC to 230 VAC

Operation:

Intermittent service: 10 mins ON / 5 mins OFF

Output characteristics:

No-load voltage: 250 Volts (nominal value without handpiece)

Min. ultrasound frequency: 28 kHz

Irrigation flow rate: from 0 to 120 ml/min (nominal value)

Adjustment in steps of 10 ml/min (nominal value)

Flush flow rate: 120 ml/min (nominal value)

Protection:

Electrical class: Class 1

Electrical safety class: BF

Safety:

Protection against internal malfunctions

2 fuses (mains connector): 5mm x 20mm / 2 AT for 100 VAC to 230 VAC

1 internal fuse that cannot be accessed by the operator reference F1: 5mm x 20mm - 10 AT / 250 VAC

Overall dimensions:

Console:

Width: 378mm

Height: 136mm

Height with bracket: 460mm

Depth: 306mm

Weight: 3700g

Pedal:

Depth: 176mm

Width: 173mm

Height with arch: 140mm

Weight: 1000g max

LCD screen:

Height: 86mm

Width: 115mm

Handpiece cable: 2000mm

Temperatures:

Operating: +10°C to +30°C

Storage: -20°C to +70°C

Humidity:

Operating: 30% to 75%

Storage: 10% to 100% including condensation

Atmospheric pressure:

Between 500 hPa and 1060 hPa

IV - INSTALLATION / COMMISSIONING

4.1 UNPACKING THE DEVICE

Upon receipt of the device, look for any damage that may have occurred during transportation.

If necessary, contact your supplier.

4.2 RECOMMENDATIONS

Ask an approved dental installation technician to make any connection to the electrical network.

The device's electrical connection shall comply with standards in force in your country.

Warning: To avoid any risk of electric shock, this device must be connected to an electrical power supply with a protective ground.

4.3 INSTALLATION

Important:

Do not install PIEZOTOME SOLO/PIEZOTOME SOLO LED near to or on another device.

Do not conceal the mains power cable and the pedal cable in a cable cover or cable guide.

- Place the console on a fixed horizontal surface or a surface that is not sloped by an angle of more than 5 degrees.
- Check that the mains isolating switch (Fig. 2-5) is in position 0 (stop).
- Connect the mains power cable to the device's mains connector (Fig. 2-1).
- Connect the mains power cable to the electrical network socket equipped with an earthing pin.
- If necessary, connect the ground equipotential plug of your installation to the potential equalizing terminal of your device (Fig. 2-6).
- Connect the control pedal to the pedal connector (Fig. 2-2).
- Position the pedal so that you can easily reach it with your foot.
- Plug the bracket into the bracket support (Fig. 2-4).

- Connect the handpiece cable to the connector (Fig. 1-7).
- Check that the device is sufficiently close to the working area to avoid having to stretch the cables. If necessary, move the device closer.
- Hang the bottle or pouch containing physiological saline solution from the bracket.

4.4 FIRST START-UP

Before using the device for the first time, it is vital that you service and/or sterilize all equipment in accordance with the procedures defined in section 11.

V - IRRIGATING SOLUTIONS

The device is not designed to administer medical substances and shall only be used with bottles or pouches containing no more than 1 litre of physiological saline solution or sterile water.

VI - IRRIGATION LINES

The device shall be used with SATELEC irrigation lines.

VII - SETTINGS / PARAMETERS / INTERFACE

7.1 START-UP

 **Important:**

It is recommended that you wait 4 seconds between stopping and restarting the device.

At start-up, the device performs an auto-check sequence and the screen displays programme d1 details.

7.2 PARAMETERS

The parameters are stored during the configuration of each programme and are recalled each time you select a programme.

7.3 FLUSH / PRIME

The flush/prime function is controlled using the keyboard.

Press button  (Fig. 3-3) to activate the flush/prime function.

This function remains active while the button is pressed down.

Flush can be activated during use of the device. The flow rate during flush is 120ml/min.

7.4 IRRIGATION

a) Setting the irrigation flow rate

Set the irrigation flow rate by pressing the buttons  and  (Fig. 3-1 and 2).

Each press adjusts the flow rate by + or - 10ml/min. The value set is then displayed on the LCD screen (Fig. 3-7) in addition to a bar graph.

The flow rate is stored in the programme being used. The flow rate can be altered during use of the device.

b) Deactivating irrigation

To deactivate irrigation, press the  button (Fig. 3-2) until 0ml/min is displayed. The irrigation bar graph (Fig. 3, item 7) is then blank.

7.5 ULTRASOUND FUNCTION

 **Important:**

The irrigation parameters can be customised by the user for each programme.

The **ULTRASOUND** function can be adjusted when the handpiece is not connected to the device.

Programme d1 is the most powerful setting. Select the programme to be used by pressing buttons

, ,  or  (Fig. 3-5) or by pressing the symbol  on the pedal (Fig. 4-2)

The selected programme appears on the LCD screen as d1, d2, d3 or d4 in addition to an associated bar graph (Fig. 3-8)

Adjust the irrigation flow rate for each programme according to the recommendations of the tip-related instructions (see chap. 7.4)

Press the  button to validate the new settings (Fig. 3-4).

The storage of parameters is confirmed by an audible and visual signal (logo  changes colour).

Repeat this sequence to set all programmes, as required.

7.6 TOOL KIT

Important:

The device has a *TOOL KIT* function used to make various adjustments, such as brightness or return to factory default settings.

To activate the *TOOL KIT* function, switch off the device, wait for 4 seconds, switch the device back on and press  and  (Fig. 3-1 and 2) as soon as they appear on the home screen (Fig. 6).

A quick press on button  (Fig. 3-4) allows you to select the brightness settings menu  (Fig.5-1) or the default factory settings menu  (Fig.5-2). The selected menu is identified by a cross (Fig. 5-4).

a) Adjusting screen brightness

The user can adjust the brightness of the device screen.

Select the brightness settings menu  (Fig.5-1).

Adjust the brightness using the buttons  and  (Fig. 3-1 and 2). The brightness can be adjusted between 0% and 100%.

The brightness intensity level is indicated (Fig. 5-5).

b) Return to default factory settings

Returning to the default factory settings solely allows irrigation flow rates to be restored.

Select the default factory settings menu  (Fig.5-2).

The various factory configurations will be effective once you exit the tool kit.

The device is configured at the factory according to the following parameters:

PROGRAMME		IRRIGATION FLOW RATES	MAINS FUNCTIONS
Very powerful	D1	60ml/min	Osteotomy, osteoplasty
Powerful	D2	60ml/min	
Average power	D3	60ml/min	
Low power	D4	60ml/min	Soft tissue removal

c) Timer

The device has a function that indicates the duration of use of the device. The total duration is displayed on screen (Fig. 5-7) under the logo  (Fig.5-3). The shortest duration displayed is one minute.

d) Software version

The software version is displayed at the bottom of the screen (Fig. 5-6).

e) Storing settings

To store the modified parameters and use the device, hold down button  for a while (Fig. 3-4). The screen shown in Fig. 6 reappears.

VIII - SAFETY

The device is fitted with a system that can detect malfunctions that may occur.

When an internal operating error occurs at start-up, symbol 1 (Chap XVII - Symbols) lights up (Fig. 7-1) and 4 beeps sound.

When an internal operating error occurs during operation, symbol 1 (Chap XVII - Symbols) lights up (Fig. 8-1) and 4 beeps sound. We recommend that you

shut down the device with the mains isolating switch (Fig. 2-5) then switch the device back on.

When the ultrasound function is activated and the handpiece is not connected or is not properly connected, symbol 2 (Chap XVII - Symbols) lights up (Fig. 9-1) and 4 beeps sound. We recommend that you connect the handpiece or check that it is correctly connected to the device.

IX - USING THE DEVICE

Important:

- Do not disconnect the handpiece cable when the device is on and the pedal is pressed down.
- Do not tighten or loosen the tips when the handpiece is activated.
- All accessories used shall be previously cleaned, disinfected and sterilized.
- For your safety and the safety of your patient, the device shall only be used with accessories supplied or recommended by SATELEC.
- Before and after each use, check the device and its accessories to make sure there are no problems.
- If necessary, do not use the device and replace any faulty parts.

Operate the device as follows:

- Make sure that the device is correctly plugged in and corrected insulated.
- Set the mains isolating switch (Fig.2-5) to position I (ON).
- Hang the bottle or pouch containing physiological saline solution or sterile water from the bracket.
- Lift the lid of the irrigating pump.
- Horizontally insert an irrigation line cassette into the purpose-designed housing.
- Lower the lid.
- Connect the irrigation line perforator to the bottle or pouch containing physiological saline solution or sterile water.
- Set the required parameters (irrigation flow rate, programme), see chapter VII.

Note:

Automatic storage allows you to save the display parameters used previously when switching between the different programmes.

This sequence is carried under controlled aseptic conditions for the safety of your patients. You may require the help of an assistant.

- Connect the handpiece cable to the device connector (Fig. 1-7).
- Screw the selected tip onto the handpiece using a wrench (according to the recommendations of the clinical manual).
- Attach the irrigation line to the cable using line clips.
- Attach the end of the irrigation line to the handpiece.
- Prime the irrigation system by pressing the flush button on the screen (Fig. 3-3)
- When the irrigation solution arrives at the tip of the tip, release the pedal.
- Use the device in compliance with good dentistry practice.

X - SWITCHING OFF THE DEVICE

At the end of the dental procedure:

- Remove the bottle or pouch containing physiological saline solution or sterile water from the bracket.
- Remove the irrigation line perforator from the bottle or pouch.
- Immerse the irrigation line perforator in a container of distilled water.
- Rinse the irrigation line and the handpiece by activating the flush function until the receptacle and the irrigation line are completely empty.
- Set the device to OFF (O).
- Disconnect the handpiece from the device.
- Remove the irrigation line clips.
- Disconnect the irrigation line from the handpiece and throw it in a bin used for the disposal of used medical items (irrigation line can only be used once). See chapter 11.4.

XI - SERVICING AND STERILIZATION

11.1 SERVICING THE DEVICE

 **Important:**

The device must be OFF during cleaning/disinfection procedures.

Do not use abrasive product, sprays or liquids to clean and disinfect the console.

Avoid the use of cleaning and disinfection products containing flammable agents.

If you do use such a product, make sure that all the product has evaporated and that no flammable agents remain on the device and its accessories before using the device again.

The device, the bracket and the multifunction pedal are not suitable for sterilization.

The console, the control pedal and the bracket shall be cleaned and disinfected (using alcohol, disinfection products, dentist disinfecting wipes, SEPTOL™ wipes) systematically after each use.

The console must be checked regularly to detect any problems.

It is important that you check the cleanliness of the ventilation inlets on the console to prevent abnormal heating.

11.2 SERVICING THE ULTRASOUND HANDPIECE

 **Important:**

The handpiece is not delivered sterile.

Refer to the user's manual delivered with the equipment.

11.3 SERVICING TIPS

 **Important:**

Tips are not delivered sterile.

Refer to the user's manual delivered with the equipment.

11.4 SERVICING IRRIGATION LINES

Sterile irrigation line

The disposable irrigation lines supplied by SATELEC shall be systematically disposed off after each use.

The reuse of a disposable irrigation line may result in contamination of your patients for which you are liable.

Do not attempt to resterilize disposable irrigation lines. Do not attempt to modify irrigation lines.

Sterilizable irrigation line

Refer to the instructions delivered with the equipment.

To service, clean, decontaminate and sterilize the sterilizable irrigation line, refer to the specific technical maintenance guide.

XII - CHECKING / MAINTENANCE

 **Important:**

In the event of an anomaly, we recommend you contact the supplier of your device. Do not contact any other repairer who may render your device hazardous for you and your patients.

12.1 CHECKING

Check the device and its accessories regularly to detect any insulation fault or any damage. Replace if necessary.

It is important to check the cleanliness of the ventilation grilles on the console ventilation to prevent abnormal overheating.

12.2 MAINTENANCE

The device does not require any preventive maintenance other than the checking of accessories and frequent cleaning/disinfection.

12.3 REPLACING FUSES

Your device is protected by two fuses in the mains connector (Fig. 2-1).

To replace these fuses, follow the steps below:

- Switch off the device (position 0).
- Disconnect the mains power cable from the electrical network.
- Disconnect the mains power cable from the mains connector (Fig. 2-1).
- Insert the tip of a flathead screwdriver into the notch located above the fuse holder to release it.
- Remove the used fuses.
- Replace the used fuses with fuses of the same type and the same amperage.
- Place the fuse holder in its recess by pushing it until you hear a click confirming it is correctly installed.
- Connect the mains power cable to the mains (Fig. 2-1).
- Connect the mains power cable to the electrical network.

Note:

The device also has an internal fuse not accessible to the user. Contact SATELEC after-sales service (see Chapter 2 - Repairs). All the information necessary for the repair of faulty components is made available by SATELEC to the technical personnel of the SATELEC-approved dealer network who carry out the work.

12.4 TROUBLESHOOTING

Refer to the table below:

Anomalies observed	Possible causes	Solutions
No operation (LCD screen off)	Faulty connection of the mains power cable	Check the mains plug. Return to the SATELEC after-sales team.
	Mains isolating switch in position 0.	Set the mains isolating switch to position I.
	No mains voltage.	Call out an electrician.
	Fuse(s) in the mains connector not working.	Replace the fuses in the mains connector.
	Internal fuse not working.	Return to the SATELEC after-sales team.
No operation (LCD screen on)	Transmission fault.	Switch off the device and switch it back it. Return to the SATELEC after-sales team.
No spray	Pouch or bottle of irrigating solution empty.	Replace the pouch or bottle of irrigating solution.
	Irrigation deactivated.	Increase the irrigation flow.
	Irrigation line blocked.	Replace the irrigation line.
	Irrigation line pinched.	Check the length of the irrigation line.
Unsuitable spray	Irrigation flow rate wrongly set.	Reset the irrigation flow rate.
Fault concerning the low vibration power of tips	Tip worn or distorted.	Replace the tip.
	Incorrect power setting.	Refer to the clinical guide for further information.
	Incorrect use: angle of approach incorrect or inadequate pressure on the clinical site.	Refer to the clinical guide for further information.
	Presence of liquid or humidity on the handpiece connector.	Thoroughly dry the electric contacts
No ultrasound	Tip not tight enough.	Retighten the tip using the spanner.
	Faulty connector contact.	Clean the connector contacts.
	Handpiece cable wire cut.	Return to the SATELEC after-sales team to replace the cable.
Liquid leak in the irrigation pump	Pipe ruptured in the irrigation line cassette.	Replace the irrigation line.
No light (Only on the PIEZOTOME SOLO LED handpiece)	LED ring not fitted to the handpiece	Fit the LED ring
	Faulty LED ring	Replace the LED ring
	Faulty LED ring connector contacts	Clean the LED ring connector contacts
	LED ring poles reversed	Position the LED ring according to the «+» foolproofing marks.
	Faulty handpiece and/or cord connector contacts	Clean the handpiece and/or cord connector contacts
	Other	Return to SATELEC after-sales

XIII - ELECTROMAGNETIC COMPATIBILITY

Important:

The power cord, the ultrasonic handpiece cords and the multifunction footswitch cord must be kept away from each other.

The device requires special precautions to be taken with regard to electromagnetic compatibility.

It must be installed and prepared for use as described in chapter 4.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the device.

The recommended separation distances in this paragraph must therefore be complied with.

The device must not be used near or on top of another device.

If this cannot be avoided, its operation under the conditions of use must be checked beforehand.

The use of accessories other than those specified or sold by SATELEC as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the device.

13. 1 ELECTROMAGNETIC EMISSIONS

The device is intended for use in the electromagnetic environment specified in the table below. The user and/ or installer must ensure that the device is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emission - CISPR 11.	Group 1	The device uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment.
RF emission - CISPR 11.	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic current emission IEC 61000-3-2.	Class A	
Voltage fluctuation and flicker IEC 61000-3-3.	Complies	

13. 2 ELECTROMAGNETIC IMMUNITY

The device is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the device is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2.	± 6 KV contact ± 8 KV air	± 6 KV contact ± 8 KV air	Floors must be wood, concrete, cement or tiled. If floors are covered with synthetic material (carpet, etc.), the relative humidity must be at least 30%.
Electrical fast transients IEC 61000-4-4.	± 2KV for power supply lines	± 2KV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment (hospital, clinic).
Surges IEC 61000-4-5.	± 1KV differential mode ± 2KV common mode	± 1KV differential mode ± 2KV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations IEC 61000-4-11.	<5% U_T (>95% dip in U_T) for 0.5 cycles. 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 250 cycles	<5% U_T (>95% dip in U_T) for 0.5 cycles. 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the use of the device requires continued operation during power mains interruption, it is recommended that the product be powered from a separate power supply (UPS, etc.).

13. 3 ELECTROMAGNETIC IMMUNITY / RADIOFREQUENCY MOBILE EQUIPMENT

The device is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the device is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Portable and mobile radiofrequency communications devices must not be used near the device (including its cables) at a distance less than that recommended and calculated according to the frequency and power of the emitter.			
Conducted disturbance, radiofrequency fields. IEC61000-4-6	3 V/m 150 KHz to 80 MHz	3 V/m	Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated radiofrequency electromagnetic field. IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80MHz to 800MHz. $d = 2.3 \sqrt{P}$ 800MHz to 2.5GHz. Where P is the maximum power rating of the emitter in watts (W) according to the manufacturer's specifications and d is the recommended minimum separation distance in meters (m).
The electromagnetic field strengths of fixed radiofrequency emitters, as determined by an electromagnetic environment measurement (a), must be less than the compliance level in each frequency range (b). Interference may occur near equipment marked with the symbol below: 			

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

(a): The electromagnetic field strengths of fixed radiofrequency emitters, such as base stations for mobile telephones (cellular/cordless), mobile radios, amateur radio, AM/FM radio broadcasts and TV broadcasts cannot be determined exactly by theory. To assess the electromagnetic environment due to fixed radiofrequency emitters, an electromagnetic environment measurement must be made. If the measured radiofrequency field strength in the immediate environment where the product is used exceeds the compliance level specified above, the performance of the product must be tested to verify whether it conforms to the specification. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.

(b): In the 150 kHz to 80 MHz frequency range, the electromagnetic field strengths must be less than 3 V/m.

13. 4 RECOMMENDED SEPARATION DISTANCES

The device is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled.

The device user and/or installer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radiofrequency communications equipment emitters and the device, according to the maximum output power of the equipment, as recommended in the table below.

Rated max. power of the emitter (W)	Separation distance in metres (m) according to emitter frequency		
	150 KHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For emitters rated at max. power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the emitter, where P is the max. power rating of the emitter in watts (W) according the manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

13. 5 CABLE LENGTHS

Cables and accessories	Maximum length	Complies with:
Handpiece cords Footswitch cord Power supply cord	Less than 3m	RF emission, CISPR 1 - Class A/Group 1 Harmonic current emission: IEC61000-3-2 Voltage fluctuation: IEC61000-3-3 Immunity to electrostatic discharge: IEC61000-4-2 Immunity to electrical fast transients/bursts: IEC61000-4-4 Immunity to surges: IEC61000-4-5 Immunity to voltage dips, short interruptions and voltage variations: IEC61000-4-11 Immunity to conducted disturbances induced by radiofrequency fields: IEC61000-4-6 Immunity to radiated radiofrequency electromagnetic fields: IEC61000-4-3

XIV - DISPOSAL AND RECYCLING

As electrical and electronic equipment, the device must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (in particular on the European market, with reference to Directive 2002/96/EC of 27/01/2003).

When your device reaches the end of its life, we consequently recommend that you contact your dental equipment dealer (or, failing this, the ACTEON Group website, the address of which is given in chapter 18), for information on how to proceed.

XV - MANUFACTURER'S LIABILITY

The manufacturer is not liable if:

- the manufacturer's installation recommendations have not been followed (supply voltage, electromagnetic environment, etc.)
- repairs or other work have been done on the device by persons not authorized by the manufacturer
- the device has been used connected to an electrical system that does not comply with current regulations
- the device has been used in ways other than those specified in this manual
- accessories (tips, handpiece, irrigation lines, etc.) other than those supplied by SATELEC have been used
- the instructions in this document have not been followed.

Note: The manufacturer reserves the right to modify the device and/or the user's manual without notice.

XVI - REGULATIONS

This medical device is classified as class IIa according to European Directive relevant to Medical Devices currently in force.

This equipment is manufactured in compliance with the current IEC 60601-1 standard.

This equipment has been designed and manufactured according to an ISO 13485-certified quality assurance system.

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