

 Via Roma, 45 21057 Olgiate Olona (VA)	Project Name: X-MIND trium	XMt Maintenance check list All. 7 a PR 7.5.3_1
		<b>Rev 1</b>




---

## PREVENTIVE MAINTENANCE CHECK LIST

---

<b>X MIND trium data</b>	
SN _____	UDI _____
Facility name _____	
Facility address _____	
State/province _____	City _____
Nation _____	Zip code _____
Facility phone number _____	E-mail _____



## PLEASE NOTE

Carry out these verifications and maintenance operations every twelve months.

After the completion of the maintenance and verification of the X-MIND trium, the installer technician MUST fill the form contained in this document to certify that the device has been correctly verified.

*The maintenance Checklist Form must be completely filled in all the details, signed and stamped by the installer technician and must be uploaded in the section "REGISTER YOUR PRODUCT" on*

<https://www.acteongroup.com>

**Upload PM Checklist, *scan or photo acceptable - must be legible***

Make a copy of the filled form and leave it to the customer (keep it with the device documentation)



## PLEASE NOTE

If you encounter problems that don't allow to correctly pass the tests, or you have any doubt for the correct maintenance of the equipment, contact immediately your referring technician or the manufacturer: ***imaging.italysupport@acteongroup.com***

<b>Test</b>	<b>Description</b>	<b>Check box</b>
<b>X-MIND TRIUM</b>		
Labels integrity	Verify that the labels are not damaged and are still visible.	<input type="checkbox"/>
Manuals	Verify the manuals contained in the specific folder on the workstation (3D versions) or in the USB key (PAN, CEPH versions) are still valid.  Please find the last version of the manuals on the web-site <a href="http://www.acteongroup.com">www.acteongroup.com</a>	<input type="checkbox"/>
U-ARM Lubrication	Check the U-Arm motions (X, Y, S) and rotation (R). Verify they are smooth, and no faults are presents.	<input type="checkbox"/>
	Lube the motors axis with DS-ES lubricant.	<input type="checkbox"/>
	Lube retaining ring and ball bearing seat of the kinematic group with OKS475 lubricant.	<input type="checkbox"/>
Tubehead visual check	Verify if oil leakage or defects present.	<input type="checkbox"/>
Mainboard fan visual check	Verify possible dust accumulation on fan of the mainboard.	<input type="checkbox"/>
COLUMN	Ensure that the bolts of wall plate are correctly tightened. Tighten if needed.	<input type="checkbox"/>
	Ensure regular column motion without jams and noises. Lube motor axis with the DS-ES lubricant.	<input type="checkbox"/>
	Check activation of mechanical limit switches of the COLUMN.	<input type="checkbox"/>
Covers visual inspection	Ensure that all the covers are intact, cover the X-MIND Trium and verify that all screws are present and fixed. Once done verify again all movements of the X-MIND Trium.	<input type="checkbox"/>
Power line values	Check power line values, measure the voltage and check that the line is dedicated to the Trium.	<input type="checkbox"/>

<b>Test</b>	<b>Description</b>	<b>Check box</b>
<b>WORKSTATION</b>		
Fans check	Verify the fans inside the workstation: if accumulated dust is present clean with compressed air if needed.	<input type="checkbox"/>
SOFTWARE	Check that no unauthorized software was installed by customer.	<input type="checkbox"/>
	Check free space on HDD (if space is limited, proceed by deleting older projections if authorized by the customer)	<input type="checkbox"/>
	Perform and save on external memory devices the backup of Trium settings and configuration with AIS Backup Manager after quality checks are passed.	<input type="checkbox"/>
<b>Test</b>	<b>Description</b>	<b>Check box</b>
<b>GENERIC MAINTENANCE OPERATIONS</b>		
LED	Ensure proper operation of the yellow LED on the control panel during x-ray emission.	<input type="checkbox"/>

EXTERNAL LAMP LIGHTS	Ensure that the external lamp lights up when the Equipment is in emission status (if lamp is present).	<input type="checkbox"/>
BUZZER	Ensure that the buzzer properly works during x-rays emission.	<input type="checkbox"/>
REMOTE EXPOSURE	Ensure that the remote exposure and emergency buttons work properly: try to release the exposure button before the acquisition is complete and ensure that the system stops, and an error message is displayed on the monitor.	<input type="checkbox"/>
LASERS	Verify patient positioning laser operation and calibration. Adjust lasers if needed.	<input type="checkbox"/>

### **QUALITY ASSESSMENT CHECKS**

*Perform the next operations with XMIND Trium completely covered with the plastics.*

*See Quality assessment manual*

<b>PAN</b> Quality Check → Absence of artifacts	Verify that the exposed area is smooth and without artifacts	<input type="checkbox"/>
<b>PAN</b> Quality Check → High Contrast Resolution	Verify that the High Contrast resolution is at least 3.1 Lp/mm	<input type="checkbox"/>
<b>PAN</b> Quality Check → Low Contrast Resolution	Verify the Low Contrast resolution, checking that all the four holes are visible.	<input type="checkbox"/>
<b>CEPH</b> Quality Check → Absence of artifacts	Verify that the exposed area is smooth and without artifacts	<input type="checkbox"/>
<b>CEPH</b> Quality Check → High Contrast Resolution	Verify that the High Contrast resolution is at least 3.1 Lp/mm	<input type="checkbox"/>
<b>CEPH</b> Quality Check → Low Contrast Resolution	Verify the Low Contrast resolution, checking that all the four holes are visible.	<input type="checkbox"/>
<b>CBCT 80x80</b> Quality Check → Homogeneity	Verify that the Homogeneity of the 80x80 FOV Test is in the acceptance range $15.0 \pm 5$ Measured value (SW output) _____	<input type="checkbox"/>
<b>CBCT 80x80</b> Quality Check → CNR	Verify that the CNR of the 80x80 FOV Test is $> 5.0$ Measured value (SW output) _____	<input type="checkbox"/>
<b>CBCT 80x80</b> Quality Check → $V_{10\%}$	Verify that the $V_{10\%}$ of the 80x80 FOV Test is $> 1.6$ LP/mm Measured value (SW output) _____	<input type="checkbox"/>
<b>CBCT 80x80</b> Quality Check → $V_{50\%}$	Verify that the $V_{50\%}$ of the 80x80 FOV Test is $> 0.5$ LP/mm Measured value (SW output) _____	<input type="checkbox"/>
<b>CBCT 110x80</b> Quality Check → Homogeneity	Verify that the Homogeneity of the 110x80 FOV Test is in the acceptance range $12.0 \pm 5$ Measured value (SW output) _____	<input type="checkbox"/>
<b>CBCT 110x80</b> Quality Check → CNR	Verify that the CNR of the 110x80 FOV Test is $> 4.0$ Measured value (SW output) _____	<input type="checkbox"/>

