



## ATTENZIONE

Utilizzare esclusivamente il lubrificante al silicone fornito da SATELEC® per lubrificare gli O-ring

## WAARSCHUWING

Gebruik uitsluitend siliconenvet van SATELEC® om o-ringen te smeren

## ΠΡΟΕΙΔΟΠΟΙΗΣΗ


Χρησιμοποιείτε αποκλειστικά το λιπαντικό σιλικόνης που παρέχεται από τη SATELEC® για τη λίπανση των δακτυλίων κυκλικής διατομής

## ATENÇÃO

Utilize apenas massa lubrificante de base de silicone fornecida pela SATELEC® para lubrificar juntas tóricas

## VARNING

Använd endast smörjfett av silikon från SATELEC® för att smörja O-ringarna



## WARNING

Exclusively use silicone grease provided by SATELEC® to lubricate o-rings

## ATTENTION


Utilisez exclusivement la graisse silicone fournie par SATELEC® pour lubrifier les joints

## WARNUNG

Für das Schmierien der Ringdichtungen ausschließlich von SATELEC® geliefertes Silikonfett verwenden

## ATENCIÓN

Utilice exclusivamente grasa de silicona proporcionada por SATELEC® para lubrificar las juntas tóricas



Use silicone grease  
Do not use mineral oil

