








RING COLOR	CLINIC	TIPS	RECOMMENDED POWER	RECOMMENDED IRRIGATION in % on NEWTRON P5XS
GREEN	Periodontics 	P2L	3	33
		P2R	3	33
		TK1-1S	2	45
		TK1-1L	2	53
		TK2-1L	2	33
		TK2-1R	2	33
		PH1	2	40
		PH2L	2	40
		PH2R	2	40
		H1	2	28
		H2L	2	28
		H2R	2	28
		H3	2	28
		H4L	2	40
		H4R	2	40
YELLOW	Endodontics 	ETBD	7	53
		ET18D	10	65
		ET20	7	48
		ET25	7	68
		ET25S	7	65
		ET25L	7	58
		ET20D	7	48
		ET40	7	73
		ET40D	7	73
		SO4	7	0
		IRR20/21	9	45
		IRR20/25	9	45
		IRR25/21	9	45
		IRR25/25	9	45
		K 15/21	6	45
		K 25/25	6	63
		CAP1	10	38
		CAP2	10	38
		CAP3	10	38
		Surgical endodontics 	AS3D	7
	AS6D		7	23
	AS9D		6	23
	ASLD		7	48
	ASRD		7	48
	P14D		7	23
	P15LD		7	48
	P15RD		7	48
	Conservatrice and restorative dentistry 	PM2	7	48
		PM3	10	50
		PMS2	7	48
PMS3		10	50	
BLUE	Prophylaxis 	1	14	43
		1S	14	43
		2	14	43
		3	14	43
		10P	14	43
		10X	12	43
		10Z	12	43
	Conservatrice and restorative dentistry 	EX1	12	43
		EX2	12	43
		EX3	12	43
		EXL	12	43
		EXR	12	43
		PM1	15	43
		PMS1	15	43
		PM4	14	43
ORANGE	Conservatrice and restorative dentistry 	5AE	20	83
		ETPR	20	83